

Pacific Health Options Inc.

Specification Documents

SCENAR NT Device Specification Parameters




World Leaders in Non-Invasive Pain Management

Drs. John Haché and Dr. Lorraine Vanbergen
Pacific Health Options Inc.

1540 10th Street East, Courtenay, B.C. V9N 7Z2 Canada

Ph: 250-338-1885 - Email: thedoctors@pacifichealthoptions.com



World Leaders in SCENAR Sales and Training

SCENAR NT DEVICE SPECIFICATIONS

Device parameters, versions	SCENAR-2-NT	SCENAR-1-NT				
		03	02 (Variant 3)	02 (Variant 2)	02 (Variant 1)	01
Diagnostics:						
EPD	-	-	-	-	-	+
IDM	-	Integral dose and null	Integral dose and null	Integral dose and null. Differential dose	Integral dose and null. Differential dose	Integral dose and null. Differential dose
Screening	-	-	-	-	1 modes	2 modes
Dose indication	Sound	Sound and symbolic	Sound and symbolic	Sound and symbolic	Sound and symbolic	Sound and symbolic
Damping	-	-	4 constant modes, 1 variable mode	4 constant modes, 1 variable mode	4 constant modes, 1 variable mode	4 constant modes, 1 variable mode
Amplitude modulation	Saw-tooth	"3 : 1"	"1 : 1", "2 : 1", "3 : 1", "4 : 1", "5 : 1"	"1 : 1", "2 : 1", "3 : 1", "4 : 1", "5 : 1"	"1 : 1", "2 : 1", "3 : 1", "4 : 1", "5 : 1"	"1 : 1", "2 : 1", "3 : 1", "4 : 1", "5 : 1", "BEE
Frequency modulation	30 - 120Hz, cycle - 8sec	30 - 120Hz, cycle - 8sec	30 - 120Hz, cycle - 8sec	30 - 120Hz, cycle - 8sec, swing	30 - 120Hz, cycle - 8sec, swing	30 - 120Hz, cycle - 8sec, swing
Influencing impulse frequencies	60, 120Hz	15, 60, 90, 350Hz	15...350Hz	15...350Hz	15...350Hz	15...350Hz
Intensity (the number of impulses in a packet)	1	1	1 ... 8	1 ... 8	1 ... 8	1 ... 8
The duration of gaps between impulses in a packet - Z	-	-	200 ... 1600µsec	200 ... 1600µsec	200 ... 1600µsec	200 ... 1600µsec
Power Influence	10 ... 250 units	1 ... 250 units	1 ... 250 units	1 ... 250 units	1 ... 250 units.	1 ... 250 units.
Indicator	Light-emitting diodes	Graphic LCD with backlight	Graphic LCD with backlight	Graphic LCD with backlight	Graphic LCD with backlight	Graphic LCD with backlight
Input of personal sets of parameters	-	-	Up to 5 sets are saved.	Up to 5 sets are saved.	Up to 5 sets are saved.	Up to 5 sets are saved.
Function of musical modulations frequency	-	Musical frequency modulations only available in SCENAR-1-NT-M	Musical frequency modulations only available in SCENAR-1-NT-M	Musical frequency modulations only available in SCENAR-1-NT-M	Musical frequency modulations only available in SCENAR-1-NT-M	Musical frequency modulations only available in SCENAR-1-NT-M